

## Signia App Remote Control Functions

The various remote control functions offered by the Signia app help you in personalizing your hearing experience further.

The main menu will show the current hearing program on the top and the current volume level of the hearing aids.

In case you wear two hearing aids, you can use either one single volume slider or two separate volume sliders, depending on the adjustment done by your Hearing Care Professional.

If you have a tinnitus program, you will find an additional icon at the bottom to adjust the volume of the tinnitus therapy signal.

By tapping on the related icon at the bottom, you can adjust the sound balance between bass and treble.

By tapping on the corresponding icon at the bottom, you can see the streaming devices in case you do audio streaming via StreamLine TV or StreamLine Mic.

FOR MORE INFORMATION: see the interactive user guide: <https://www.signia-hearing.com/signia-app/>



If you are using an audio streaming accessory, e.g. StreamLine TV, you will be able to adjust its volume level via this menu.

The Directional Hearing function in the Signia app offers you the possibility to adjust the directionality of your hearing aids according to your needs.

By tapping the icon on the right side at the bottom, you can activate the Directional Hearing functionality.



Directional Hearing is normally set to automatic. The "Auto" button at the bottom part is marked red in this case.

You can adjust the Directional Hearing manually by selecting its span. In this case, the selected span is marked red while the "Auto" button turns white.

Directional Hearing can be adjusted flexibly on the 360° basis depending on your listening environment and conversation partner.

← → ← → ← → ← →

Directional Hearing can be adjusted flexibly on the 360° basis depending on your listening environment and conversation partner.

Directional Hearing can be adjusted flexibly on the 360° basis depending on your listening environment and conversation partner.

The Signia app enables you to monitor the connection and battery status of your hearing aids in an intuitive way.

← → ← → ← → ← →

The Bluetooth and battery icons show your connection status and battery level. By tapping on one of these icons, you will get a more detailed information.

You can monitor intuitively if the Bluetooth connection between the app and the hearing aids is lost. In order to check the battery status, you can tap on one of the battery icons shown on the menu.

A sound played by your hearing aids will indicate the battery status:  
1 beep = low  
2 beeps = medium  
3 beeps = full

For assistance with using the Signia app, please call the **Consumer Hotline: 800-350-6093**